


Astrovirus Strip

One Step Astrovirus Antigen Test Strip

A rapid, one step test for the qualitative detection of *Astrovirus* antigens in human feces.

IVD For In-Vitro diagnostic and professional use only



INTENDED USE

The *Astrovirus* Strip is a rapid chromatographic immunoassay for the qualitative detection of *Astrovirus* antigens in human faeces specimens to aid in the diagnosis of *Astrovirus* infection.

INTRODUCTION

Viral gastroenteritis is an infection caused by a variety of viruses that results in vomiting or diarrhea. Many different viruses can cause gastroenteritis, including rotaviruses, noroviruses, adenoviruses, sapoviruses, and astroviruses.

The main symptoms of viral gastroenteritis are watery diarrhea and vomiting. The affected person may also have headache, fever, and abdominal cramps ("stomach ache"). In general, the symptoms begin 1 to 2 days following infection with a virus that causes gastroenteritis and may last for 1 to 10 days, depending on which virus causes the illness. Some research studies have shown that the duration of the symptoms are approximately three to four days. Astrovirus infection is not usually a severe situation and only in some rare cases leads to dehydration. Adenoviruses and astroviruses cause diarrhea mostly in young children, but older children and adults can also be affected.

PRINCIPLE

The *Astrovirus* Strip is a qualitative immunoassay for the detection of *Astrovirus* antigen in human faeces samples. The membrane is pre-coated with monoclonal antibodies against *Astrovirus* antigens on the test line region. During testing, the sample reacts with the particle coated with anti-*Astrovirus* antibodies which was pre-dried on the test strip. The mixture moves upward on the membrane by capillary action. In the case of a positive result the specific antibodies present on the membrane will react with the mixture conjugate and generate a coloured line. A green coloured band always appears in the control line and serves as verification that sufficient volume was

added, that proper flow was obtained and as an internal control for the reagents.

PRECAUTIONS

- For professional in vitro diagnostic use only.
- Do not use after expiration date.
- The test should remain in the sealed pouch until use.
- Do not use the test if pouch is damaged.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, do not eat, drink or smoke in the area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The test should be discarded in a proper biohazard container after testing.
- The test must be carried out within 2 hours of opening the sealed bag.

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at refrigerated or room temperature (2-30°C/36-86°F). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

MATERIALS

MATERIALS PROVIDED

- Test Strips
- Buffer
- Package Insert

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer
- Disposable gloves

SPECIMEN COLLECTION AND PREPARATION

Collect sufficient quantity of faeces (1-2 g or mL for liquid sample). Stool samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator (2-4°C/36-40°F) for 1-2 days prior to testing. For longer storage the specimen must be kept frozen at -20°C/-4°F. In this case, the sample will be totally thawed, and brought to room temperature before testing.

PROCEDURES

To process the collected stool samples (see illustration 1):

Use a separate specimen collection vial/testing tube for each sample. Dispense exactly 1mL of the buffer into a vial. Introduce the swab or stick two times into the faecal specimen to pick up quite a lot of sample (200-300mg). and put into the testing tube

or vial with buffer. Shake the testing tube or vial in order to assure good sample dispersion. For liquid stool samples, aspirate the faecal specimen with a dropper and add 200-300 µL into the testing tube or vial with buffer.

Test Procedure (see illustration 2)

Allow the tests, stool samples and buffer to reach to room temperature (15-30°C/59-86°F) prior to testing. Do not open pouches until ready to perform the assay.

1. Remove the *Astrovirus* Strip from its sealed pouch and use it as soon as possible.
2. Shake the specimen collection vial to assure a good sample dispersion. Break off the cap of the vial.
3. Use a separate test strip for each sample. Leave the test strip to stand vertically taking care of not surpassing the limit of immersion indicated with the arrows. Start the timer.
4. Read the result at 10 minutes after dispensing the sample.

Illustration 1






Illustration 2



INTERPRETATION OF RESULTS

Illustration 3



POSITIVE: Two lines appears across the central window in the result line region (**red** test line marked in the illustration 3 with

the letter T) and in the control line region (**green** control line marked in the illustration 3 with the letter C).

NEGATIVE: Only one **green** band appears across the control line region marked in the illustration 3 with the letter C (control line).

INVALID: A total absence of the green control coloured band regardless the appearance or not of the red test line. Note: Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (T) will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test:

- A green line appearing in the control line region (C). It confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. *Astrovirus* Strip will only indicate the presence of *Astrovirus* in the specimen (qualitative detection) and should be used for the detection of *Astrovirus* antigens in faeces specimens only. Neither the quantitative value nor the rate of increase in *Astrovirus* antigens concentration can be determined by this test.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
3. Some stool samples can decrease the intensity of the control line.
4. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of *Astrovirus* infection.
5. This test provides a presumptive diagnosis of *Astrovirus* infections. All results must be interpreted together with other clinical information and laboratory findings available to the physician.

EXPECTED VALUES

A study of intestinal disease in the UK, published in 1999 determined incidence as 3.8/1000 patient years in the community (95%CI, range 2.3-6.4), the fourth most common known cause of viral gastroenteritis. Studies in the USA have detected astroviruses in the stools of 2-9% of children presenting symptoms; illness is most frequent in children of less than two years, although outbreaks among adults and the elderly have been reported. Early studies carried out in Glasgow demonstrated that a significant proportion of babies excreting virus particles, 12%, did not exhibit gastrointestinal symptoms, and seroprevalence studies carried out in the US have shown that 90% of children have antibody to *HastV-1* by age 9, suggesting that (largely asymptomatic) infection is common. There is, as with most viral causes of gastroenteritis, a peak of incidence in the winter.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

The evaluation was conducted comparing the results obtained using the *Astrovirus* Strips to a commercial available *Astrovirus* ELISA assay.

The detection of *Astrovirus* showed >94% of concordance in sensitivity and >99% of concordance in specificity.

Cross-Reactivity

It was performed an evaluation to determine the cross reactivity of *Astrovirus* Strip. There is not cross reactivity with common gastrointestinal pathogens, other organisms and substances occasionally present in faeces.

















- Rotavirus
- Adenovirus
- *Escherichia coli*
- *Campylobacter*
- *Giardia lamblia*
- Human Hemoglobin

REFERENCES

1. SUNITA SHASTRI *et al.*, "Prevalence of Astroviruses in a Children's Hospital", *JOURNAL OF CLINICAL MICROBIOLOGY* Sept. 1998, p. 2571-2574 Vol. 36, No. 9.
2. ASHLEY *et al.*, "Astrovirus-associated gastroenteritis in children", *Journal of Clinical Pathology*, 1978, 31, 939-943.

 **Atlas Medical**
Ludwig-Erhard Ring 3
15827 Blankenfelde-Mahlow
Germany
Tel: +49 - 33708 – 3550 30
Email: Info@atlas-medical.com

PPI1838A01
Rev A (02.09.2019)

	Catalogue Number		Temperature limit
	In Vitro diagnostic medical device		Caution
	Contains sufficient for <n> tests and Relative size		Consult instructions for use (IFU)
	Batch code		Manufacturer
	Do not re-use		Use-by date
	Manufacturer fax number		Do not use if package is damaged
	Manufacturer telephone number		Date of Manufacture
	Keep away from sunlight		Keep dry