

## Atlas One Step Cocaine Test Strip (Urine)

A rapid, one step test for the qualitative detection of Cocaine metabolites in human urine.

**IVD** For *in vitro* diagnostic use only.

Store at (2-30°C)

### INTENDED USE

Atlas One Step Cocaine Test Strip (Urine) is a rapid chromatographic immunoassay for the qualitative detection of Cocaine metabolite, Benzoyllecgonine, in human urine at a cut-off concentration of 300 ng/mL.

This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

### INTRODUCTION

Cocaine, is a potent central nervous system (CNS) stimulant and a local anesthetic. Initially, it brings about extreme energy and restlessness while gradually resulting in tremors, over-sensitivity and spasms. In large amounts, Cocaine causes fever, unresponsiveness, and difficulty in breathing and unconsciousness.

Cocaine is often self-administered by nasal inhalation, intravenous injection and free-base smoking. It is excreted the urine in a short time primarily as Benzoyllecgonine. Benzoyllecgonine, a major metabolite of Cocaine, has a longer biological half-life (5 - 8 hours) than Cocaine (0.5 - 1.5 hours), and can generally be detected for 24-48 hours after Cocaine exposure.

Atlas One Step Cocaine Test Strip (Urine) is a rapid urine-screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Cocaine metabolite in urine. Atlas One Step Cocaine Test Strip (Urine) yields a positive result when the Cocaine metabolite in urine exceeds 300 ng/mL. This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA).

### PRINCIPLE

Atlas One Step Cocaine Test Strip (Urine) is an immunoassay based on the principle of competitive binding.

Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody.

During testing, a urine specimen migrates upward by capillary action. Benzoyllecgonine, if present in the urine specimen below 300 ng/mL, will not saturate the binding sites of antibody in the test strip. The antibody-coated particles will then be captured by immobilized Benzoyllecgonine conjugate and a visible colored line will appear in the test line region. The colored line will not form in the test line region if the Benzoyllecgonine level is above 300 ng/mL because it will saturate all the binding sites of antibodies.

A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

### MATERIALS

#### Materials Provided

- Test strips (Contain mouse monoclonal anti-Benzoyllecgonine antibody-coupled particles and Benzoyllecgonine-protein conjugate. A goat antibody is employed in the control line system).
- Package insert

#### Materials Required But Not Provided

- Specimen collection container
- Timer
- External controls

### PRECAUTIONS

- For professional *in vitro* diagnostic use only. Do not use after the expiration date.
- The test strip should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test strip should be discarded in a proper biohazard container after testing.

### STORAGE AND STABILITY

- The kit can be stored at room temperature or refrigerated (2-30°C).
- The test Cassette is stable through the expiration date printed on the label on the sealed pouch.
- The test Cassette must remain in the sealed pouch until use. Do not freeze.
- Do not use beyond the expiration date.

### SPECIMEN COLLECTION AND PREPARATION

#### Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens

exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.


### Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

### PROCEDURE

**Allow the test strip, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.**

1. Bring the pouch to room temperature before opening it. Remove the test strip from the sealed pouch and use it as soon as possible.
2. With arrows pointing toward the urine specimen, **immerse the test strip vertically in the urine specimen for at least 10-15 seconds.** Do not pass the maximum line (MAX) on the test strip when immersing it. See the illustration below.
3. Place the test strip on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear. The result should be **read at 5 minutes.** Do not interpret the result after 10 minutes.



**INTERPRETATION OF RESULTS** (Please refer to the illustration above)

**NEGATIVE: Two lines appear.** One red line should be in the control region (C), and another apparent red or pink line should be in the test region (T). This negative result indicates that the Benzoyllecgonine concentration is below the detectable level (300 ng/mL).

**NOTE:** The shade of red in the test region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

**POSITIVE: One red line appears in the control region (C).** No line appears in the test region. This positive result indicates that the Benzoyllecgonine concentration is above the detectable level (300 ng/mL).

**INVALID: Control line fails to appear.** Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a

new test strip. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

### QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

### LIMITATION

1. Atlas One Step Cocaine Test Strip (Urine) provides only a qualitative, preliminary analytical result. A secondary quantitative analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.
2. It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
3. Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
4. A positive result does not indicate level of intoxication, administration route or concentration in urine.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. Test does not distinguish between drugs of abuse and certain medications.

### PERFORMANCE CHARACTERISTICS

#### Accuracy

A side-by-side comparison was conducted by laboratory personnel using Atlas One Step Cocaine Test Strip (Urine) and a commercially available rapid test. Testing was performed on specimens previously collected from subjects presenting for Drug Screen Testing. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method		Other COC Rapid Test		Total Results
Results	Positive	Negative		
Atlas One Step Cocaine Test Strip	Positive	136	0	136
	Negative	7	157	164

Total Results	143	157	300
% Agreement with this commercial kit	95%	>99%	98%

When compared to GC/MS at the cut-off of 300 ng/mL, the following results were tabulated:

Method		GC/MS		Total Results
Results	Positive	Negative		
Atlas One Step Cocaine Test Strip	Positive	119	17	136
	Negative	5	159	164
Total Results		124	176	300
%Agreement with GC/MS Analysis		96%	90%	93%

#### Analytical Sensitivity

A drug-free urine pool was spiked with Benzoylcegonine at the following concentrations: 0 ng/mL, 150 ng/mL, 225 ng/mL, 300 ng/mL, 375 ng/mL and 450 ng/mL. The result demonstrates 100% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

Benzoylcegonine Concentration (ng/mL)	Percent of Cut-off	n	Visual	Result
			Negative	Positive
0	0	30	30	0
150	-50%	30	30	0
225	-25%	30	30	0
300	Cut-off	30	4	26
375	+25%	30	0	30
450	+50%	30	0	30

#### Analytical Specificity

The following table lists compounds that are positively detected in urine by Atlas One Step Cocaine Test Strip (Urine) at 5 minutes.

Compound	Concentration (ng/ml)
Benzoylcegonine	300
Cocaine HCl	780
Cocaethylene	12500
Ecgonine HCl	32000

#### BIBLIOGRAPHY

1. Stewart DI, T Inoba, M Ducassen, W Kalow. Clin. Pharmacol. Ther. 1979; 25:264
2. Ambre J. J. Anal. Toxicol. 1985; 9:241
3. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man.

2nd Ed. Biomedical Publ., Davis, CA. 1982: 488 Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986



**ATLAS MEDICAL**  
**Ludwig-Erhard Ring 3**  
**15827 Blankenfelde-Mahlow**  
**Germany**  
**Tel: +49 - 33708 – 3550 30**  
**Email: Info@atlas-medical.com**  
**PPI1693A01**  
**Rev A (02.09.2019)**

	Catalogue Number		Temperature limit
	In Vitro diagnostic medical device		Caution
	Contains sufficient for <n> tests and Relative size		Consult instructions for use (IFU)
	Batch code		Manufacturer
	Do not re-use		Use-by date
	Manufacturer fax number		Do not use if package is damaged
	Manufacturer telephone number		Date of Manufacture
	Keep away from sunlight		Keep dry